



ASPEN IMAGING
HEALTHCARE

ASPEN VetTable



Digital X-Ray
for Veterinary

ASPEN Vet



ASPENVet





Table



SYSTEM SPECIFICATIONS

STAND	4-way Floating Stand, Fixed Table
TUBE (E-7239X)	Target Angle 16 , 125kVp, 140kHU, Focal Spot (1.0 / 2.0)
COLLIMATOR	MC-D (M)
GENERATOR (POSKOM, 32kW)	High Frequency type Touch Screen Control mA Range: 10 - 400 mA kVp Range: 40 - 125 kVp mAs Range: 0.1-500 mAs Single Phase 220 V
HIGH VOLTAGE CABLE	40 VA
MAXIMUM WEIGHT CAPACITY	Up to 150 Kg

DIMENSIONS (mm)





ASPEN IMAGING
HEALTHCARE

801 Jupiter Rd Suite 200, Plano, TX 75074, USA
sales@aspenimaging.com | (+1) 214-257-0113

WWW.ASPEN-IMAGING.COM

DISTRIBUTED BY

About Aspen Imaging

Aspen Imaging Healthcare is a Texas-based Original Equipment Manufacturer (OEM) providing complete imaging solutions for x-ray distributors and custom designed, private-labeled solutions for OEMs across the Americas.

Founded in 2015, Aspen Imaging has become one of the fastest growing Medical Imaging companies with offices in the United States, Mexico and South Korea.

