



ASPEN IMAGING
HEALTHCARE

Europa85L


Portable X-ray System Unit



CLINICALLY ENGINEERED FOR HIGH MOBILITY AND EFFICIENT IMAGING WORKFLOWS—SUPPORTING RADIOLOGY TEAMS WHEREVER CARE IS NEEDED

The Europa85L is Aspen Imaging's compact yet powerful portable X-ray solution, engineered for mobility, precision, and everyday reliability. With an ultra-light weight of just 8 lb, it's ideal for mobile healthcare teams and small practice environments that require fast setup, effortless transport, and dependable performance. Designed to deliver consistent high-quality imaging without compromising workflow efficiency, the Europa85L supports clinicians who need flexibility, durability, and premium image results in any setting—from bedside exams to field diagnostics.

DEVICE FEATURES



* When using Mirror View, the image transmitted via Miracast is not intended for diagnostic purposes.

**Available under specific conditions.

WHAT'S INCLUDED

- Europa85L
- Carrying case
- Handswitch
- Charger
- Carrying strap



EXTREME LIGHT WEIGHT

8 lb for ultimate
portability and convenience

SPECIFICATIONS Europa85L

POWER OUTPUT	425 W
kV / mA	50 ~ 85 kV / 2~5 mA
EXPOSURE TIME	0.03 ~ 2 sec
COLLIMATOR	Multi-leaf collimator
ANODE TYPE	Stationary anode
ANODE ANGLE	16°
FOCAL SPOT	0.8 mm
HEAT CAPACITY	8,000 J
BATTERY CAPACITY	22.2 VDC / 2,200 mAh
OPERATION	LCD touch screen
WEIGHT	About 3.6 kg / 8 lb
DIMENSIONS	266 mm (W) x 167 mm (H) x 156 mm (D)

HANDHELD MODE

