



ASPEN IMAGING
HEALTHCARE

laica90

Portable X-ray System Unit from the AspenVet line.



laica90

What's included

- Hard Case
- Hand Switch
- Bluetooth Mouse
- Touch pen



TECHNICAL SPECIFICATIONS

POWER OUTPUT	450 W
kV / mA	50 ~ 90 kV / 2 ~ 5 mA
EXPOSURE TIME	0.04 ~1.3 sec
COLLIMATOR	Multi-leaf collimator
ANODE TYPE	Stationary anode
TARGET ANODE ANGLE	19°
FOCAL SPOT	0.5 mm
ANODE HEAT STORAGE	10,000 J
BATTERY CAPACITY	22.2 VDC / 2,200 mAh
WEIGHT	4.9 kg ± 5% / 9.9 lb
DIMENSIONS	367 x 201 x 188 mm / 14.40" x 7.87" x 6.31"



LAICA PC SPECIFICATIONS

SYSTEM INFORMATION	
MACHINE NAME	Laica
OPERATING SYSTEM	Windows 11 Pro 64 bits
CPU	Intel Gemini Lake J4125 @ 2.00 GHz (4 CPUs), ~2.0 GHz
RAM	8 GB
SSD	128 GB
DISPLAY SPECS	
CARD NAME	Intel(R) HD Graphics
MANUFACTURER	Intel Corporation
SCREEN SIZE	10.1" LCD touch screen
SCREEN RESOLUTION	1920X1200 (max)
PORTS	
WIFI MODULES	2(2.5-5 GHz)
BLUETOOTH MODULE	1
USB PORTS	2: 1(2.0), 1 (3.0)

COMPATIBLE RADIOGRAPHIC SYSTEM



LAICA STAND STANDARD



ASPEN IMAGING
HEALTHCARE

801 Jupiter Rd Suite 200, Plano, TX 75074, USA
sales@aspenimaging.com | (+1) 214-257-0113

WWW.ASPEN-IMAGING.COM

DISTRIBUIDO POR

About Aspen Imaging

Aspen Imaging Healthcare is a Texas-based Original Equipment Manufacturer (OEM) providing complete imaging solutions for x-ray distributors and custom designed, private-labeled solutions for OEMs across the Americas.

Founded in 2015, Aspen Imaging has become one of the fastest growing Medical Imaging companies with offices in the United States, Mexico and South Korea.

