

# T4 OTC Ceiling X-Ray System

PRODUCT BOOKLET



ASPEN IMAGING  
HEALTHCARE



# A Work of Engineering Art

The T4 OTC is our ceiling system that elegantly combines a sophisticated, feature-rich rad room with an intuitive and instinctive operation for your technologists. Developed with functionality and convenience in mind, we've meticulously thought through every detail of our T4 to ensure the system is a pleasure to work with for years to come.

- ✓ HOSPITAL GRADE QUALITY
- ✓ ERGONOMICALLY DESIGNED
- ✓ FLEXIBLE SYSTEM



# Device Features

## A.I. TUBE CRANE

Sturdy 5-tier telescoping column that allows for head to toe examination possibilities

## ERGONOMIC HANDLE

Ergonomic handles with “all lock release button” for soft handling

## ADVANCED MANUAL COLLIMATOR

Manual collimator with LED light (adjustable brightness) and a bright, red-laser crosshair

## PILLAR-STYLE ELEVATING TABLE

Pillar-style elevating table on both sides

## MOTORIZED WALL STAND

Motorized wall stand uses electromagnetic braking for smooth movement

## OPTIONAL HANDLES

Wall stand is accompanied by optional overhead handle and patient side handle

## MEDICAL GRADE TABLETOP

Medical grade white tabletop with optional transparent tabletop (in picture) that allows the operator to see detector alignment with the tube







# Tube Crane Display

Position your system toward the table or wall stand and have the correct receptor highlighted.

Ergonomic operator handles with an “all lock release” button that allows the operator to smoothly glide the tube crane with minimal physical effort.

Integrated with live feedback of SID and degrees.

## Stylish, Yet Rugged Elevating Table

Ultra-high-end elevating table using actuators (not motors) for a smooth and silent operation.


Standard tabletop provides a white surface to match the rest of the system – optional transparent tabletop provides a modern design, while allowing the operator to be able to ensure that the detector is always aligned to the tube. Includes a wired foot switch that can be moved to different areas, depending on the preference of the operator.





# Integration with Viewworks Detectors & Software

- Powerful image processing engine with simple user interface.
- Anytime™: Viewworks' proprietary AED (Automatic Exposure Detection) technology
- Track post-exam statistics with dose and exam management tools
- Rugged: 100 cm drop tested flat panel detectors Medical Imaging companies with offices in the United States, Canada, Mexico and South Korea.



**VIEWWORKS**

## About Aspen Imaging

Aspen Imaging Healthcare is a Texas-based Original Equipment Manufacturer (OEM) providing complete imaging solutions for x-ray distributors and custom designed, private-labeled solutions for OEMs across the Americas.

Founded in 2015, Aspen Imaging has become one of the fastest growing Medical Imaging companies with offices in the United States, Canada, Mexico and South Korea.





ASPEN IMAGING  
HEALTHCARE

801 Jupiter Rd Suite 200, Plano, TX 75074, USA  
sales@aspensate.com | (+1) 214-257-0113

[WWW.ASPEN-IMAGING.COM](http://WWW.ASPEN-IMAGING.COM)